
Medtrum EasyView Webové stránky

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

 English

[Contact Us](#) | [Privacy Policy](#) | [Terms of Use](#)
Copyright © 2023 Medtrum Technologies Inc. All rights reserved.

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English

Medtrum

Apply for an HCP account

United Kingdom (mmol/l)

Name

Organization

Enter email address

Enter mobile number

I Agree to the [Terms of Use of Medtrum](#)

I Agree to the [Privacy Policy of Medtrum](#)

I Agree NOT to share the data to third parties without patient's consent

Submit

Back

Contact Us | Privacy Policy | Terms of Use
Copyright © 2022 Medtrum Technologies Inc. All rights reserved.

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
 English

[Contact Us](#) | [Privacy Policy](#) | [Terms of Use](#)
 Copyright © 2023 Medtrum Technologies Inc. All rights reserved.

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
The screenshot shows the 'Patients' list view in the Medtrum application. It includes a 'Delete' button and a 'Find Patient' search bar. The table below lists patient records with columns for Name, Gender, Age, Type, Upload time, and Action.

ID	Name	Gender	Age	Type	Upload time	Action
ID	Jack	Male		Type 1		
ID	Thomas					
ID	Maria	Male				
ID	Linda	Male	40	Type 1		
ID	Peterson					
ID	James					
ID	Robert	Male		Type 1		
ID	John					

Showing 1 to 8 of 15 entries

Footer: Contact Us | Privacy Policy | Terms of Use | Copyright © 2022 Medtrum Technologies Inc. All rights reserved.

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


<input type="checkbox"/>	Name	Gender	Age	Type	Upload time	Action
<input checked="" type="checkbox"/>	Jack	Male		Type 1		 
<input checked="" type="checkbox"/>	Thomas					 
<input type="checkbox"/>	Maria	Male				 
<input type="checkbox"/>	Linda	Male	40	Type 1		 
<input type="checkbox"/>	Peterson					 
<input type="checkbox"/>	James					 
<input type="checkbox"/>	Robert	Male		Type 1		 
<input type="checkbox"/>	John					 

Showing 1 to 8 of 15 entries

Contact Us | Privacy Policy | Terms of Use

Copyright © 2022 Medtrum Technologies Inc. All rights reserved.

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
Delete Find Patient

<input type="checkbox"/>	Name	Gender	Age	Type	Upload time	Action
<input type="checkbox"/>	Demo Data					✎ ✖
<input type="checkbox"/>	Demo12	Male				✎ ✖

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
Data To View

Uploaded Data
The data is from EasyTouch App or EasyView iPad App.

Real-time Data
The data is from EasySense App or EasyPatch App.

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		Glucose Target Range		Time Period	Post Meal Analysis
		Low	High	hours	hours
Show alerts		OFF		Show Glucose Value OFF	
Very High		13.9		High 10.0	
Low		3.9		Very Low 3.0	
Hypoglycaemia		3.0			
Breakfast	Before	3.9	7.2	07:00 - 09:00	01:00 - 03:00
	After	5.0	8.8		
Lunch	Before	3.9	7.2	12:00 - 13:00	01:00 - 03:00
	After	5.0	8.8		
Dinner	Before	3.9	7.2	18:00 - 19:00	01:00 - 03:00
	After	5.0	8.8		
Night-time		5.0	8.8	23:00 - 06:00	

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
BG		CGM		Insulin		Carbs		Exercise	
Average		Average		Average Daily Dose		Average Carbs/Day		Average time/Day	
9.5 mmol/L		9.6 mmol/L		53.58 U		5 g		-- min	
SD	#	SD	#	SD	# days	SD	#	SD	#
5.48	18	2.61	4680	8.10	3.0	0.0	1	--	--


Estimated A1C: 7.7% calculated from SG values Note: Estimated A1C does not replace lab measurement.

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
BG	BG / calibration (per day)	4.5 / 2.5
	AUC high > 5.6 mmol/L	3.8
	AUC low < 3.9 mmol/L	0
	Highest value	22.2 mmol/L
	Lowest value	2.2 mmol/L
SG	<3.1 mmol/L	0.1%
	<3.9 mmol/L	0.3%
	3.9 mmol/L-5.6 mmol/L	2.7%
	>5.6 mmol/L	97%
	>13.9 mmol/L	12.1%
Sensor	Sensor used	3
Pump	Daily Basal	48.63U/90.8 %
	Daily Bolus	4.95U/9.2 %
	Average Bolus No. (per day)	0.3
	Average Suspend time (per day)	18 min
	Average Reservoir usage (per day)	0.72 days


Manual Bolus (Volume/No. per day)	0U/0
-Normal	0U/0
-Extended	0U/0
-Combo	0U/0
Bolus Calculator (Volume/No. per day)	0.07U/0.3
-Food Bolus calc.	0.07U/0.3
-Corr Bolus calc.	0U/0
-Override(+)	0U/0
-Override(-)	0U/0

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Auto Mode	Time in Auto Mode	86%
	Average SG ± SD	8.3 ± 2.77 mmol/L
	Time in Target Range	73%
	Time above Range	25%
	Time below Range	2%
	Auto Basal (per day)	12.89U
Auto Bolus (per day)	11.24U	
Auto Mode Exits	No Sensor readings	0
	Max Delivery	0
	Min Delivery	0

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
Dashboard **Daily Summary** Day by Day Overview Sensor Overlay Trend Analysis Events Logbook Alert

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
Statistics						
	Glucose		Insulin		Events	
	BG	SG				
Avg. (mmol/L)±SD	7.2±1.33	6.6±1.19	Total daily dose (U)	16.95	Carbs (g)	40
Min (mmol/L)	5.6	4.3	Bolus (U)/%	1.10/6%	# of meals	2
Max (mmol/L)	8.5	9.3	Basal (U)/%	15.85/94%	Avg±SD	20.0±0.0
Target (%)	100	100	Food Bolus (U)	0	Total exercise time	30
Low (%)	0	0	Correction Bolus (U)	0	# of exercises	1
High (%)	0	0	Food+corr (U)	0.65	Avg. Intensity	30
# of readings	4	418	Normal Bolus (U)	0.65		
# of hypo/duration(min)	0	0/0	Extended Bolus (U)	0		
# of alerts	--	2	Temp Basal (U)	0		
			Duration of suspension(min)	0		
			# of pump alarms	0		

Auto Mode	
Time in Auto Mode(%)	23%
Avg.(mmol/L)±SD	7.1 ± 1.16 mmol/L
Target Range(%)	100%
Above Range(%)	0%
Below Range(%)	0%
# of Auto Mode Exits	0

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Dashboard Daily Summary **Day by Day Overview** Sensor Overlay Trend Analysis Events Logbook Alert

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
The screenshot shows the 'Create' dialog box overlaid on the application interface. The dialog contains the following fields and controls:

- Time:** Two input fields, one for 'Date' and one for 'Time'.
- Event Type:** A dropdown menu with 'Select' as the current value.
- Details:** A text input field.
- Notes:** A text input field.
- Source:** A dropdown menu with 'Web' selected.
- Buttons:** A green 'Create' button and a red 'Cancel' button.

The background interface shows a navigation menu with options like 'Dashboard', 'Daily Summary', 'Day by Day Overview', 'Sensor Overlay', 'Trend Analysis', 'Events', 'Logbook', and 'Alert'. The 'Events' section is active, showing a list of events with columns for Time, Event Type, Details, Notes, Source, and Action.

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


The screenshot shows a form for selecting the duration and date of the report. It includes a 'Duration' field with a dropdown menu set to '7 days' and a 'Date' field with a date range from '03-12-2022' to '08-12-2022'. Below the date field is a 'Week' field with a dropdown menu set to 'Mon'.

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4.5.8 Výstrahy

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Dashboard Daily Summary Day by Day Overview Sensor Overlay Trend Analysis Events Logbook **Alert**

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Duration: **7 days** 14 days 30 days 90 days Customize

Type: **Pump** CGM PDM

Date: 01-03-2023 - 07-03-2023

Time	Type	Alert
07-03-2023 21:52:10	Pump	Predictive Low Suspend
07-03-2023 19:14:10	Pump	Predictive Low Suspend
06-03-2023 21:59:00	Pump	Auto Mode Exit - No Readings

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Dashboard Daily Summary Day by Day Overview Sensor Overlay Trend Analysis Events Logbook Alert

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Bolus		Maximum Bolus: 25U		Bolus Calculator: On		Active Insulin Time: 3H	
IC Ratios(g/U)							
Start	00:00	07:30	11:00	12:30	16:30	19:30	
Carbs	7.8	6.7	7.2	8	7.5	7	
Insulin Sensitivity mmol/L							
	Start					00:00	
	BG					2.8	
BG Target							
	From					00:00	
	Range					5.0-6.1	
Bolus Reminder (Off)							
	From	--	--	--	--	--	
	To	--	--	--	--	--	

Pump			Patch Expiration: Off			Low Reservoir: Insulin 12U		
Pump Base SN: 3F2D4GT0			Command Reminder: Off			Auto Off: On 12H		
Out of Range: Off			Hourly Max: 35 U					
Daily Max: 80U								
Alarms & Alerts on APP			Alarms on App			Alerts on App		
Audio	Vibrate	Pump Light	Repeat	Snooze	Tone	Repeat	Snooze	Tone
On	Off	On	5	60	Default	10	60	Default

CGM			Sensor Expiration: Off			Glucose Alerts: Off		
Transmitter SN: A1B2C3D4								
Factory Calibration: On								
Glucose Limits								
	From					00:00		
	Low					4.2		
	High					10.0		

Predictive Alerts (h/hmm)			Rate Alerts (mmol/L/min)					
High	On	30	Rise	On	0.220			
Low	On	30	Fall	On	0.220			
Alerts for Low			Alerts for High			Alerts for Others		
Repeat	2		Repeat	2		Repeat	5	
Snooze	15		Snooze	30		Snooze	30	
Tone	Default		Tone	Notification		Tone	Notification	
Alerts for Hypo								
Tone	Notification							

APP		
Name: Bob Papgeorgiou	Audio Options: Audio	Glucose Unit: mmol/L
Weight Unit: kg	Height Unit: cm	
EasyLoop		
Low Suspend: Off	Predictive Low Suspend: Off	
Auto Mode		
Target SG: 5.6 mmol/L	Auto Meal Handling: On	Weight: 65 kg
Total Daily Dose: 62 U		Auto Mode: Off

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

Profile:			
Username: abc044	Name: abc044	Type of diabetes: Unspecified	
Height: 175 cm	Weight:	Age:	
Duration:			
Complications:			Create
Concomitant diseases	Date of diagnosis	Suggestions	Action
No Data			
Prescription:			Create
Medicine	Dose	Notes	Action
No Data			
Other suggestions:			Create
06-12-2022 - 12-12-2022			
HCP	Date	Suggestions	Action
No Data			

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Select duration: **One Week** One Month One Quarter

End date: **Yesterday** Last Saturday Last Sunday Customize

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

Number	User	Time	Operation Type	Patient
1	Jason	09-12-2022 14:43:58	Add	abcd44
2	Mariam	06-07-2022 08:43:16	Delete	Demo 12

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

Zprávy, profil a heslo

Security Settings

Multi-factor authentication (MFA) off
Every time you log in, in addition to entering a password, you also need to enter a dynamic password to further enhance the security of your account.

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